


**Product Name: Recombinant Human PYDC5 Protein**  
**Catalog #: PEH10198**



## Summary

|                              |                                                                                                                                                                                                                 |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | PYDC5_HUMAN                                                                                                                                                                                                     |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                             |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                                |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                           |
| <b>Accession #</b>           | W6CW81                                                                                                                                                                                                          |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                          |
| <b>Species</b>               | Human                                                                                                                                                                                                           |
| <b>Molecular Mass</b>        | 11 kDa                                                                                                                                                                                                          |
| <b>Buffer</b>                | Lyophilized from Phosphate buffered saline                                                                                                                                                                      |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                              |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | LTSLDNITDE ELDRFKCFLP DEFNIATGKL HTLNSTSSQL DLKRWHGVCS<br>EEDRIFQKLN YMLVAKCLRE EQETGICGSP                                                                                                                      |

## image



## Background

**Alternative Names**      Pyrin domain-containing protein 5 (Pyrin domain-only protein 3)

## Note

For Research Use Only , Not for Diagnostic Use