


Summary

Name	RGPD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	PODJDO
Host	E.coli
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ISEPGNQEKK SEKPLENDTG FQAQDISGQK NGRGVIFGQT SSTFTFADVA KSTSGEGFQF GKKDPNFKGF SGAGEKLFSS QCGKMANKAN TSGDFEKDDD ACKTEDSDDI HFEPVVQMPE KVELVTGEEG EKVLYSQGVK LFRFDAEISQ WKERGLGNLK ILKNEVNGKP RMLMRRDQVL KVCANHWITT TMNLKPLSGS DRAWMWLASD FSDGDAKLER LAAQFKTPEL AEEFKQKFEE CQRLLLDIPL QTPHKLVDTG RAAKLIQRAE EMKSGLKDFK TFLTNDQTKV TEEENKSGST GAAGASDTTI KPNPENTGPT

image



Background

Alternative Names	RANBP2-like and GRIP domain-containing protein 1 (Ran-binding protein 2-like 6) (RanBP2-like 6) (RanBP2L6)
--------------------------	--

Product Name: Recombinant Human RGD1 Protein
Catalog #: PEH10201



Note

For Research Use Only , Not for Diagnostic Use