


Product Name: Recombinant Human SCRN2 Protein
Catalog #: PEH10204



Summary

Name	SCRN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q96FV2
Host	E.coli
Species	Human
Molecular Mass	67 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YIEVEQVSKT HAVILSRPSW LWGAEMGANE HGVCIGNEAV WTKEPVGEGE ALLGMDLLRL ALERSSSAQE ALHVITGLLE HYQGQGNCLE DAAPFSYHST FLLADRTEAW VLETAGRLWA AQRIQEGARN ISNQLSIGTD ISAQHPELRT HAQAKGWWDG QGAFDFAQIF SLTQQPVRME AAKARFQAGR ELLRQRQGGI TAEVMMGILR DKESGICMDS

image



Background

Alternative Names Secernin-2

Note

Product Name: Recombinant Human SCRIN2 Protein
Catalog #: PEH10204



For Research Use Only , Not for Diagnostic Use