


Product Name: Recombinant Human SCRN2 Protein
Catalog #: PEH10206



Summary

Name	SCRN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q96FV2
Host	E.coli
Species	Human
Molecular Mass	86 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KNSDRPRDEV QEVVFVPAGT HTPGSRLQCT YIEVEQVSKT HAVILSRPSW LWGAEMGANE HGVCIGNEAV WTKEPVGEGE ALLGMDLLRL ALERSSSAQE ALHVITGLLE HYGQGGNCLE DAAPFSYHST FLLADRTEAW VLETAGRLWA AQRIQEGARN ISNQLSIGTD ISAQHPELRT HAQAKGWWDG QGAFDFAQIF SLTQQPVRME AAKARFQAGR ELLRQRQGGI TAEVMMGILR DKESGICMDS GGFRTTASMV SVLPQDPTQP CVHFLTATPD PSRSVFKPFI FGMGVAQAPQ VLSPTFGAQD PVRTLPRFQT QVDRRHLYR GHQAALGLME RDQDRGQQLQ QKQQDLEQEG LEATQGLLAG EWAPPLWELG SLFQAFVKRE

image



Background

Alternative Names	Secernin-2
--------------------------	------------

Product Name: Recombinant Human SCRIN2 Protein
Catalog #: PEH10206



Note

For Research Use Only , Not for Diagnostic Use