


Product Name: Recombinant Human SCRN2 Protein
Catalog #: PEH10207



Summary

Name	SCRN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q96FV2
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QEVVFVPAGT HTPGSRLQCT YIEVEQVSKT HAVILSRPSW LWGAEMGANE HGVCIGNEAV WTKEPVGEGE ALLGMDLLRL ALERSSSAQE ALHVITGLLE HYGQGGNCLE DAAPFSYHST FLLADRTEAW VLETAGRLWA AQRIQEGARN ISNQLSIGTD ISAQHPELRT HAQAKGWWDG QGAFDFAQIF SLTQQPVRME AAKARFQAGR ELLRQRQGGI TAEVMMGILR DKESGICMDS

image



Background

Alternative Names Secernin-2

Note

Product Name: Recombinant Human SCRIN2 Protein
Catalog #: PEH10207



For Research Use Only , Not for Diagnostic Use