


Product Name: Recombinant Human SAMD15 Protein
Catalog #: PEH10213



Summary

Name	SAM15_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P1V8
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKILELPDET KPRETHVEFS KEDRPEPIKS KYSVGNDLE HREPKRGKLS LSDKFRKEYY ALGSLRESEE SIGTHYEFLQ PLQKLLNVSE ECSYSDPSES QTEELSEFVHE KEVVDLSQEL KERVSEDDT QPEKGTQLQF EHLNWDPEEV AEWISQLGFP QYKECFITNF ISGRKLIHVN CSNLPQMGIT NFEDMKAISR HTQELLEIEE PLFKRSISLP YRDIIGLYLE QKGHTGIKSD SLTLSEFVKA AGLQDYAPEI

image



Background

Alternative Names	Sterile alpha motif domain-containing protein 15 (SAM domain-containing protein 15)
--------------------------	---

Note

Product Name: Recombinant Human SAMD15 Protein
Catalog #: PEH10213



For Research Use Only , Not for Diagnostic Use