


Product Name: Recombinant Human TTC31 Protein
Catalog #: PEH10217



Summary

Name	TTC31_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q49AM3
Host	E.coli
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TGQSDRGKGA EGLGTYCGLR KSFLYPPQES EPCPQSPSAS ATFPSVSDSL LQVAMPQKLL VTEEEANRLA EELVAEEERM KQKAEEKRLK KKRQKERKRQ ERLEQYCGEP KASTTSDGDE SPPSSPGNPV

image



Background

Alternative Names Tetratricopeptide repeat protein 31 (TPR repeat protein 31)

Note

For Research Use Only , Not for Diagnostic Use