


Product Name: Recombinant Human METTL2A Protein
Catalog #: PEH10219



Summary

Name	MET2A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96IZ6
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ERKVVQENSIQ RVCQEKQVDY EINAHKYWND FYKIHENGFF KDRHWLFTEF PELAPSQNQN HLKDWFLENK SEVPECRNNE DGPGLIMEEQ HKCSSKSLEH KTQTLPVEEN VTQKISDLEI

image



Background

Alternative Names tRNA N(3)-methylcytidine methyltransferase METTL2A (EC 2.1.1.-)
(Methyltransferase-like protein 2A)

Note

For Research Use Only , Not for Diagnostic Use