


Summary

Name	UBP12_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75317
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QFPVNEHYFG LVNFGNTCYC NSVLQALYFC RPFREKVLAY KSQPRKESL LTCLADLFHS IATQKKKVG V IPPKKFITRL RKENELFDNY MQQDAHEFLN YLLNTIADIL QEERKQEKQN GRLPNGNIDN ENNNSTPDPT WVHEIFQGTL TNETRCLTCE TISSKDEDFL DLSVDVEQNT SITHCLRGFS NTETLCSEYK YYCEECSKQ EAHKRMKVKK LPMILALHLK RFKYMDQLHR YTKLSYRVVF PLELRLFNTS GDATNPDRMY DLVAVVVHCG SGPNRGHYIA IVKSHDFWLL FDDDIVEKID AQAIEEFYGL TSDISKNSSES GYILFYQSRD

image



Background

Alternative Names	Ubiquitin carboxyl-terminal hydrolase 12 (EC 3.4.19.12) (Deubiquitinating enzyme 12) (Ubiquitin thioesterase 12) (Ubiquitin-hydrolyzing enzyme 1)
--------------------------	---

Product Name: Recombinant Human USP12 Protein
Catalog #: PEH10222



(Ubiquitin-specific-processing protease 12)

Note

For Research Use Only , Not for Diagnostic Use