


Product Name: Recombinant Human C12orf50 Protein
Catalog #: PEH10225



Summary

Name	CL050_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NA57
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQEGEIPQLS QSQEPLKPQE NISRPIHHPL VLKTNFEEEE EVDEQNDASS LWTKTPEEIE EKRAIKEMCY KSGEYYRFHT PPDILSSKSM TPTAEKQLEK PLENGSELQE GDSLTVPTKL SQYERQGEIK TSLHGKPKTD IAAFENGSGD CYVPQRVIFL GVDESEALTE EKEITISKCS NTKDNKDSPH PKHSLTTRLV PTTHVLNATE

image



Background

Alternative Names Uncharacterized protein C12orf50

Note

Product Name: Recombinant Human C12orf50 Protein
Catalog #: PEH10225



For Research Use Only , Not for Diagnostic Use