


Product Name: Recombinant Human C19orf47 Protein
Catalog #: PEH10229



Summary

Name	CS047_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N9M1
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QDMCKAATES VPCSPSPLAG EIRRG TSAAS RMITNSLNHD SPPSTPPRRP DTSTSKISVT VSNKMAAKSA KATAALARRE EESLAVPAKR RRVTAEMEGK YVINMPKGTT PRTRKILEQQ QAAKGLHRTS VFDR LG AETK ADTTTGSKPT GVFSRLGATP ETDEDLAWDS DNDSSSSVLQ

image



Background

Alternative Names Uncharacterized protein C19orf47

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C19orf47 Protein
Catalog #: PEH10229

