


Product Name: Recombinant Human C2orf73 Protein
Catalog #: PEH10236



Summary

Name	CB073_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N5S3
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HKIEDAAITY VSENEEIKHE EKP G KSIHHS KSHVGRGRIY YAKFINTNAR TYNEPFPYID PKKGPEIQGD WWSHGKALEP VFLPPYDSKS TQRSD FQKPS CPLVLPVKHS KMQKPSCGIV PLASPGTSAE LQNNFIEYIS FIHQYDARKT PNEPLQGKRH GAFVQREIKP GSRPTVPKGA EVLLNTPGSR SSEQSKKTEK GNSAESRMIS PGLCQQNSQE LLEPKTHLSE TDVRQAAKAC PSTPESREKT SGATQTTVGD ALFTRHKPLN

image



Background

Alternative Names Uncharacterized protein C2orf73

Product Name: Recombinant Human C2orf73 Protein
Catalog #: PEH10236



Note

For Research Use Only , Not for Diagnostic Use