


Product Name: Recombinant Human C3orf62 Protein
Catalog #: PEH10238



Summary

Name	CC062_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZUJ4
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PAFATNHAPV NAKPHALCPE RKPLTSKENV LMHSSILAPE RESWRTAGEG ENWRKENLRK DMERDLKADS NMPLNNSQEQ VTKDLLDMID HTSIRTIEEL AGKIEFENEL NHMCGHCQDS PFKEEAWALL MDKSPQKATD ADPGSLKQAF DDHNIVETVL DLEEDYNVMT

image



Background

Alternative Names Uncharacterized protein C3orf62

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C3orf62 Protein
Catalog #: PEH10238

