


Product Name: Recombinant Human C6orf132 Protein
Catalog #: PEH10239



Summary

Name	CF132_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5T0Z8
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ENTASSAPVD WRDPSQMEKL RNELAAYLCG SRREDRFLSH RPGPTVAPQS KEGKKGPRLP EKETLLSLPA KDTPPGVPEK SLGGSSLTET EAAPSLTLPS VDYIPQDSPT PSVRQIRNEL EARLSSAAEK EAKPSIGSLP PKPRLEGGRI CENGADDDKL SKPVAKNLPP

image



Background

Alternative Names Uncharacterized protein C6orf132

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C6orf132 Protein
Catalog #: PEH10239

