


Product Name: Recombinant Human C9orf131 Protein
Catalog #: PEH10242



Summary

Name	C1131_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5VYM1
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SHCHHFQSSS QLQPQESGRA EQGSQRGEKM KGKMVSQVPS QGPCVHMEAG VDYLSPGPGE PSNSKVLVSG KRKDKASASS SAKKREHPRK PKAGDHRRTG ARLGLSTVTG KNHPAQARSL VEAPVSTFPQ RSQHRGQSSQ HTALPQLLLP KASGPQDQPE

image



Background

Alternative Names Uncharacterized protein C9orf131

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human C9orf131 Protein
Catalog #: PEH10242

