


Product Name: Recombinant Human YC006 Protein
Catalog #: PEH10246



Summary

Name	YC006_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZUG5
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KDLSQEFGMP LSQEELTDEK LFALPPQPAP NLEDYHSRNS TLTLEIHAHQ EPRKRFTYSQ DYLSAMVEPL DLKEEEKKAQ KKSQRQAWLTA RGFQVTGLQS DTESSFQDLK LPPIKELNEE WKENSLFANV LEPVLDRDRW SWDRHHVDFD LYKKPPPFLE LLPSPAPKPV

image



Background

Alternative Names Uncharacterized protein FLJ43738

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human YC006 Protein
Catalog #: PEH10246

