


Product Name: Recombinant Human RIMOC1 Protein
Catalog #: PEH10250



Summary

Name	CE051_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NDU8
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RVEELGDLAQ AHQQQLSEAA GEDDHFLIRA SAALEKLKLL CGEEKECSNP SNLLELYTQA ILDMTYFEEN KLVDEDFPED SSSQKVKELI SFLSEPEILV KENNMHPKHC NLLGDELLEC LSWRRGALLY MYCHSLTKRR EWLLRKSSLL KKYLLDGISY LLQMLNYRCP IQLNEGVSFQ DLDTAKLLSA GIFSDIHLLA MMYSGEMCYW GSKYCADQQP ENHEVDTSVS GAGCTTYKEP LDFREVGEKI LKKYVSVCEG PLKEQEWNTT

image



Background

Alternative Names UPF0600 protein C5orf51

Product Name: Recombinant Human RIMOC1 Protein
Catalog #: PEH10250



Note

For Research Use Only , Not for Diagnostic Use