



Product Name: Recombinant Human MYT1L Protein
Catalog #: PEH10256



Summary

Name	MYT1L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UL68
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TVDSLKLLAQ GHGVVLSNM NDRNYADSMS QQDSRNMNYV MLGKPMNNGL MEKMVEESDE EVCLSSLECL RNQCFLARK LSETNPQERN PQQNMNIRQH VRPEEDFPGR TPDRNYSMDL NLMRLEEQLS PRSRVFASCA KEDGCHERDD DTTSVNSDRS EEVFDMTKGN LTLLEKAIAL ETERAKAMRE KMAMEAGRDR

image



Background

Alternative Names Myelin transcription factor 1-like protein (MyT1-L) (MyT1L)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MYT1L Protein
Catalog #: PEH10256

