


Product Name: Recombinant Human CPSF3 Protein
Catalog #: PEH10260



Summary

Name	CPSF3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKF6
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIQEKPAKLV FKNITVIQEP GMVVLEWLAN PSNDMYADTV TTVILEVQSN PKIRKGAVQK VSKKLEMHVY SKRLEIMLQD IFGEDCVSVK DDSILSVTVD GKTANLNLET RTVECEESE DDESLREMVE

image



Background

Alternative Names	Cleavage and polyadenylation specificity factor subunit 3 (EC 3.1.27.-) (Cleavage and polyadenylation specificity factor 73 kDa subunit) (CPSF 73 kDa subunit) (mRNA 3'-end-processing endonuclease CPSF-73)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CPSF3 Protein
Catalog #: PEH10260

