


Product Name: Recombinant Human PCF11 Protein
Catalog #: PEH10261



Summary

Name	PCF11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O94913
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSPSISTPPI VPDIQKNLTQ EQLIRQQLLA KQKQLLELQQ KKLELELEQA KAQLAVLSLV QQETSNLGPG SAPSKLHVSQ IPPMAVKAPH QVPVQSEKSR PGPSLQIQDL KGTNRDPRLN RISQHSHGKD QSHRKEFLMN TLNQSDTKTS KTIPSEKLNS SKQEKSXSGE KITKELDQL DSKSKSKSKS PSPLKNKLSH

image



Background

Alternative Names Pre-mRNA cleavage complex 2 protein Pcf11 (Pre-mRNA cleavage complex II protein Pcf11)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PCF11 Protein
Catalog #: PEH10261

