


Product Name: Recombinant Human ZMYND11 Protein
Catalog #: PEH10265



Summary

Name	ZMY11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15326
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NQSVEPKKEE PEPETEAVSS SQEPTMPQP IEKVSSTQT KKLSASSPRM LHRSTQTTND GVCQSMCHDK YTKIFNDFKD RMKSDHKRET ERVREALEK LRSEMEEEKR QAVNKAVANM QGEMDRKCKQ VKEKCKEEFV EEIKKLATQH KQLISQTKKK QWCYNCEEEA

image



Background

Alternative Names	Zinc finger MYND domain-containing protein 11 (Adenovirus 5 E1A-binding protein) (Bone morphogenetic protein receptor-associated molecule 1) (Protein BS69)
--------------------------	---

Note



Product Name: Recombinant Human ZMYND11 Protein
Catalog #: PEH10265

For Research Use Only , Not for Diagnostic Use