


Product Name: Recombinant Human NSMCE2 Protein
Catalog #: PEH10268



Summary

Name	NSE2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96MF7
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ASSVALDLVE SQTEVSSEYS MDKAMVEFAT LDRQLNHVVK AVQSTINHVK EERPEKIPDL KLLVEKKFLA LQSKNSDADF QNNEKFVQFK QQLKELKKQC GLQADREADG TEGVDEDIIV TQSQTNFTCP ITKEEMKKPV KNKVCGHTYE EDAIVRMIES RQKRKKKAYC PQIGCSHTDI RKSDLIQDEA

image



Background

Alternative Names	E3 SUMO-protein ligase NSE2 (EC 2.3.2.-) (E3 SUMO-protein transferase NSE2) (MMS21 homolog) (hMMS21) (Non-structural maintenance of chromosomes element 2 homolog) (Non-SMC element 2 homolog)
--------------------------	--

Note

Product Name: Recombinant Human NSMCE2 Protein
Catalog #: PEH10268



For Research Use Only , Not for Diagnostic Use