


Product Name: Recombinant Human ZRSR2 Protein
Catalog #: PEH10269



Summary

Name	U2AFM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15696
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RQELARLRDS GLSQKEEED TFIEEQLEE EKLLERERQR LHEEWLLREQ KAQEEFRIKK EKEEAARKRQ EEQERKLKEQ WEEQQRKERE EEEQKRQEKK EKEEALQKML DQAENELENG TTWQNPEPPV DFRVMEKDRA NCPFYSKTGA

image



Background

Alternative Names U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 2 (CCCH type zinc finger, RNA-binding motif and serine/arginine rich protein 2) (Renal carcinoma antigen NY-REN-20) (U2(RNU2) small nuclear RNA auxiliary factor 1-like 2) (U2AF35-related protein) (URP)

Note

Product Name: Recombinant Human ZRSR2 Protein
Catalog #: PEH10269



For Research Use Only , Not for Diagnostic Use