


Product Name: Recombinant Human YTHDC1 Protein
Catalog #: PEH10270



Summary

Name	YTDC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96MU7
Host	E.coli
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKKGSKRKSD RMESTDTKRQ KPSVHSRQLV SKPLSSSVSN NKRIVSTKGK SATEYKNEEY QRSERNKRLD ADRKIRLSSS ASREPYKNQP EKTCVRKRDP ERRAKSPTPD GSERIGLEVD RRASRSSQSS KEEVNSEEYG SDHETGSSGS SDEQGNNTEN EEFGVEEDVE EDEEVEEDAE EDEEVDEDGE EEEEEEEEE EEEEEEEEEEY EQDERDQKEE GNDYDTRSEA SDSGSESVSF TDGSVRSRSGS TDGSDEKKKE RKRARGISPI VFDRSGSSAS ESYAGSEKKH

image



Background

Alternative Names	YTH domain-containing protein 1 (Splicing factor YT521) (YT521-B)
--------------------------	---

Product Name: Recombinant Human YTHDC1 Protein
Catalog #: PEH10270



Note

For Research Use Only , Not for Diagnostic Use