

**Product Name: Recombinant Human GTF2F2 Protein**  
**Catalog #: PEH10271**



## Summary

<b>Name</b>	T2FB_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P13984
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	12 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	QTLTVFTESS SDKLSLEGIV VQRAECPAA SENYMLKRL QIESSKPVR LSQQLDKVVT TNYKPVANHQ YNIEYERKKK EDGKRARADK

## image



## Background

<b>Alternative Names</b>	General transcription factor IIF subunit 2 (EC 3.6.4.12) (ATP-dependent helicase GTF2F2) (General transcription factor IIF 30 kDa subunit) (Transcription initiation factor IIF subunit beta) (TFIIF-beta) (Transcription initiation factor RAP30)
--------------------------	--

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human GTF2F2 Protein**  
**Catalog #: PEH10271**

---

