



Product Name: Recombinant Human UBR2 Protein
Catalog #: PEH10273



Summary

Name	UBR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8I WV8
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VAETEGTIME ESSRDKDKAE RKRKAEIARL RREKIMAQMS EMQRHFIDEN KELFQQTLEL DASTSAVLDH

image



Background

Alternative Names	E3 ubiquitin-protein ligase UBR2 (EC 2.3.2.27) (N-recognin-2) (RING-type E3 ubiquitin transferase UBR2) (Ubiquitin-protein ligase E3-alpha-2) (Ubiquitin-protein ligase E3-alpha-II)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use