


Product Name: Recombinant Human STAG3 Protein
Catalog #: PEH10274



Summary

Name	STAG3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UJ98
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GDSLLAEDT DFEDSLNRNV KKRAAKRPPK TTPVAKHPKK GSRVVHRHSR KQSEPPANDL FNAVKAAKSD MQSLVDEWLD SYKQDQDAGF LELVNFFIQS CGCKGIVTPE MFKKMSNSEI IQHLTEQFNE

image



Background

Alternative Names Cohesin subunit SA-3 (SCC3 homolog 3) (Stromal antigen 3) (Stromalin-3)

Note

For Research Use Only , Not for Diagnostic Use