


Product Name: Recombinant Human PATZ1 Protein
Catalog #: PEH10275



Summary

Name	PATZ1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HBE1
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGPSNFCSIC NRGFSSASYL KVHVKTHHGV PLPQVSRHQE PILNGGAAFH CARTYGNKEG QKCSHQDPIE SSDSYGDLSL ASDLKTPEKQ SANGSFSCDM AVPKNKMESD

image



Background

Alternative Names POZ-, AT hook-, and zinc finger-containing protein 1 (BTB/POZ domain zinc finger transcription factor) (Protein kinase A RI subunit alpha-associated protein) (Zinc finger and BTB domain-containing protein 19) (Zinc finger protein 278) (Zinc finger sarcoma gene protein)

Note

Product Name: Recombinant Human PATZ1 Protein
Catalog #: PEH10275



For Research Use Only , Not for Diagnostic Use