


Product Name: Recombinant Human ZBTB18 Protein
Catalog #: PEH10277



Summary

Name	ZBT18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99592
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SLSQRSVTSV RDSADVDCVL DLSVKSSLG VENLNSSYFS SQDVLRSLV QVKVEKEASC DESDVGTNDY DMEHSTVKES VSTNNRVQYE PAHLAPLRED SVLRELDRED KASDDEMMTP

image



Background

Alternative Names	Zinc finger and BTB domain-containing protein 18 (58 kDa repressor protein) (Transcriptional repressor RP58) (Translin-associated zinc finger protein 1) (TAZ-1) (Zinc finger protein 238) (Zinc finger protein C2H2-171)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZBTB18 Protein
Catalog #: PEH10277

