


Product Name: Recombinant Human MTREX Protein
Catalog #: PEH10279



Summary

Name	MTREX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P42285
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGTKKDKKED KGKWKGPPTS ADKAGKRFDG KLQSESTNNG KTKRDVDFEG TDEPIFGKKP RIEESITDL SLADLMPRVK VQSVETVEGC

image



Background

Alternative Names	Exosome RNA helicase MTR4 (EC 3.6.4.13) (ATP-dependent RNA helicase DOB1) (ATP-dependent RNA helicase SKIV2L2) (Superkiller viralicidic activity 2-like 2) (TRAMP-like complex helicase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use