


Product Name: Recombinant Human BAZ1A Protein
Catalog #: PEH10280



Summary

Name	BAZ1A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NRL2
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GADVTSANAK YRYQKRGGFD ATDDACMELR LSNPSLVKKL SSTSVDLTP GEKMKILHAL CGKLLTLVST RDFIEDYVDI LRQAKQEFRE LKAEQHRKER EAAAARIRKR KEEKLKEQEQ KMKEKQEKLK EDEQRNSTAD ISIGEEERED FDTSIESKDT EQKELDQDMV TEDEDDPGSH KRGRRGKRGQ NGFKEFTRQE QINCVTREPL TADEEEALKQ

image



Background

Alternative Names	Bromodomain adjacent to zinc finger domain protein 1A (ATP-dependent chromatin-remodeling protein) (ATP-utilizing chromatin assembly and remodeling factor 1) (hACF1) (CHRAC subunit ACF1) (Williams syndrome transcription factor-related chromatin-remodeling factor 180) (WCRF180)
--------------------------	---

Product Name: Recombinant Human BAZ1A Protein
Catalog #: PEH10280



(hWALp1)

Note

For Research Use Only , Not for Diagnostic Use