



Product Name: Recombinant Human KHDRBS2 Protein
Catalog #: PEH10281



Summary

Name	KHDR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5VWX1
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAEKDSLDPF FVHASRLAE EIEKFQGSDF KKEDEEKYL DVISNKNIKL SERVLIPVKQ YPKFNFVGKL LGPRGNSLKR LQEETGAKMS ILKGSMRDK AKEEELRKSG EAKYAHLSDE LHVLIIEVFAP PGEAYSRMSH ALEEIKKFLV PDYNDEIRQE QLRELSYLNQ SEDSGRGRGI RGRGIRIAPT APSRGRGGAI PPPPPPGRGV LTPRGSTVTR GALPVPPVAR GVPTPRARGA PTPVPGYRAPP PPAHEAYEEY GYDDGYGGEY DDQTYETYDN SYATQTQSVP EYYDYGHGVS EDAYDSYAPE

image



Background

Alternative Names	KH domain-containing, RNA-binding, signal transduction-associated protein 2 (Sam68-like mammalian protein 1) (SLM-1) (hSLM-1)
--------------------------	---

Product Name: Recombinant Human KHDRBS2 Protein
Catalog #: PEH10281



Note

For Research Use Only , Not for Diagnostic Use