

Product Name: Recombinant Human MIXL1 Protein
Catalog #: PEH10286



Summary

Name	MIXL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2W2
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SRRQSGKSFQ PLARPEIILN HCAPGTETKC LKPQLPLEVD VNCLPEPNGV GGGISDSSSQ GQNFETCSPL SEDIGSKLDS

image



Background

Alternative Names	Homeobox protein MIXL1 (Homeodomain protein MIX) (hMix) (MIX1 homeobox-like protein 1) (Mix.1 homeobox-like protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use