


Product Name: Recombinant Human KDM8 Protein
Catalog #: PEH10287



Summary

Name	KDM8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8N371
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKRPARGSLP EQPCTKKARA DHGLIPDVKL EKTVPRLHRP SLQHFREQFL VPGRPVILKG VADHWPCMQR WSLEYIQEIA GCRTVPVEVG SRYTDEEWSQ

image



Background

Alternative Names	Bifunctional peptidase and arginyl-hydroxylase JMJD5 (EC 1.14.11.-) (EC 3.4.-.-) (JmjC domain-containing protein 5) (Jumonji C domain-containing protein 5) (L-arginine (3R)-hydroxylase KDM8)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use