


Product Name: Recombinant Human IRX4 Protein
Catalog #: PEH10288



Summary

Name	IRX4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P78413
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EEEAREEPLK SSKNAEPVGK EEKELESDL DDFDPLEAEP PACELKPPFH SLDGGGLRVP AAPDGPVKEA SGALRMSLAA GGGAALDEDL ERARSCLRSA AAGPEPLPGA EGGPQVCEAK LGFVPAGASA GLEAKPRIWS LAHTATAAAA AATSLSQTEF PSCMLKRQGP AAPAAVSSAP ATSPSVALPH SGALDRHQDS

image



Background

Alternative Names	Iroquois-class homeodomain protein IRX-4 (Homeodomain protein IRXA3) (Iroquois homeobox protein 4)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IRX4 Protein
Catalog #: PEH10288

