


Product Name: Recombinant Human PSMC3IP Protein
Catalog #: PEH10290



Summary

Name	HOP2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2W1
Host	E.coli
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VSDADLQVLD GKIVALTAKV QSLQQSCRYM EAELKELSSA LTTPEMQKEI QELKKECAGY RERLKNIAA TNHVTPEEKE QVYRERQKYC KEWRKRKRMA TELSDAILEG YPKSKKQFFE

image



Background

Alternative Names	Homologous-pairing protein 2 homolog (Nuclear receptor coactivator GT198) (PSMC3-interacting protein) (Proteasome 26S ATPase subunit 3-interacting protein) (Tat-binding protein 1-interacting protein) (TBP-1-interacting protein)
--------------------------	---

Note

Product Name: Recombinant Human PSMC3IP Protein
Catalog #: PEH10290



For Research Use Only , Not for Diagnostic Use