


Product Name: Recombinant Human ZRANB3 Protein
Catalog #: PEH10291



Summary

Name	ZRAB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5FWF4
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FTHFEKEKQH DIRSFFVPQP KKRQLMTSCD ESKRFREENT VVSSDPTKTA ARDIIDYESD VEPETKRLKL AASEDHCSPTS EETPSQSKQI RTPLVESVQE AKAQLTTPAF PVEGWQCCLC TYINNSLPY CEMCETPQGS AVMQIDSLNH IQDKNEKDDS QKDTSKKVQT ISDCEKQALA QSEPGQLADS KEETPKIEKE DGLTSQPGNE

image



Background

Alternative Names	DNA annealing helicase and endonuclease ZRANB3 (Annealing helicase 2) (AH2) (Zinc finger Ran-binding domain-containing protein 3) [Includes: DNA annealing helicase ZRANB3 (EC 3.6.4.-); Endonuclease ZRANB3 (EC 3.1.-.-)]
--------------------------	--

Product Name: Recombinant Human ZRANB3 Protein
Catalog #: PEH10291



Note

For Research Use Only , Not for Diagnostic Use