


Product Name: Recombinant Human AFF4 Protein
Catalog #: PEH10293



Summary

Name	AFF4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UHB7
Host	E.coli
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GYKKEGREQG TGNSYTDTSG PKETSSATPG RDSKTIQKGS ESGRGRQKSP AQSDSTTQRR TVGKKQPKKA EKAAAEPRG GLKIESETPV DLASSMPSSR HKAATKGSRK PNIKKESKSS PRPTAEKKKY KSTSKSSQKS REIIETDTSS SDSDESESLP PSSQTPKYPE SNRTPVKPSS VEEEDSFFRQ RMFSPMEEKE LLSPLSEPDD RYPLIVKIDL NLLTRIPGKP YKETEPKGE KKNVPEKHTR EAQKQASEKV SNKGKRKHKN EDDNRASESK KPKTEDKNSA GHKPSSNRES SKQSAAKEKD LLPSPAGVP SKDPKTEHGS RKRTISQSSS LKSSNSNKE TSGSSKNSSS TSKQKKTEGK TSSSSKEVKE KAPSSSSNCP PSAPTLDSSK PRRTKLVFDD RNYSADHYLQ EAKKLKHNAD ALSDRFEKAV YYLDAVVSFI ECGNALEKNA QESKSPFPMY SETVDLIKYT MKLKNYLAPD

image



Background

Product Name: Recombinant Human AFF4 Protein
Catalog #: PEH10293



Alternative Names

AF4/FMR2 family member 4 (ALL1-fused gene from chromosome 5q31 protein) (Protein AF-5q31) (Major CDK9 elongation factor-associated protein)

Note

For Research Use Only , Not for Diagnostic Use