


Product Name: Recombinant Human SSBP3 Protein
Catalog #: PEH10298



Summary

Name	SSBP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BWW4
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EHSSSEAKAFH DYSAAAAPSP VLGNIIPPNDG MPGGPIPPGF FQGPPGSQPS PHAQPPPHNP SSMMGPHSQP FMSPRYAGGP RPIRMGNQP PGGVPGTQPL LPNSMDPTRQ QGHPNMGGSM QRMNPPRGMG PMGPGPQNYG SGMRPPPNL GPAMPGINMG PGAGRPWPNP NSANSIPYSS SSPGTYVGPP GGGGPPGTPI MPSPADSTNS SDNIYTMINP VPPGGSRSNF PMGPGSDGPM GGMGGMEPHH MNGSLGSGDI DGLPKNSPNN ISGISNPPGT PRDDGELGGN FLHSFQNDNY

image



Background

Alternative Names	Single-stranded DNA-binding protein 3 (Sequence-specific single-stranded-DNA-binding protein)
--------------------------	---

Product Name: Recombinant Human SSBP3 Protein
Catalog #: PEH10298



Note

For Research Use Only , Not for Diagnostic Use