


Product Name: Recombinant Human L3MBTL1 Protein
Catalog #: PEH10299



Summary

Name	LMBL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y468
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NPPEDPNQDP PEDDSTCQCQ ACGPHQAAGP DLGSSNDGCP QLFQERSVIV ENSSGSTSAS ELLKPMKKRK RREYQSPSEE ESEPEAMEKQ EEGKDPEGQP TASTPESEEW SSSQPATGEK

image



Background

Alternative Names	Lethal(3)malignant brain tumor-like protein 1 (H-l(3)mbt) (H-l(3)mbt protein) (L(3)mbt-like) (L(3)mbt protein homolog) (L3MBTL1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use