



Product Name: Recombinant Human FIGLA Protein
Catalog #: PEH10303



Summary

Name	FIGLA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6QHK4
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QLVLERRVA NAKERERIKN LNRGFARLKA LVPFLPQSRK PSKVDILKGA TEYIQVLS DL LEGAKDSKKQ DPDEQSYSNN

image



Background

Alternative Names	Factor in the germline alpha (FIGalpha) (Class C basic helix-loop-helix protein 8) (bHLHc8) (Folliculogenesis-specific basic helix-loop-helix protein) (Transcription factor FIGa)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use