


Product Name: Recombinant Human DBX1 Protein
Catalog #: PEH10305



Summary

Name	DBX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NMT0
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GMLRRAVFSD VQRKALEKMF QKQKYISKPD RKKLAAKLGL KDSQVKIWFQ NRRMKWRNSK ERELLSSGGC REQTLPTKLN PHPDLSDVGQ KGPGNEEEEE GPGSPSHRLA YHASSDPQHL RDPRLPGPLP PSPAHSSSPG KPSDFSDSEE

image



Background

Alternative Names Homeobox protein DBX1 (Developing brain homeobox protein 1)

Note

For Research Use Only , Not for Diagnostic Use