


Product Name: Recombinant Human RRP1 Protein
Catalog #: PEH10306



Summary

Name	RRP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P56182
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRLRLDKFYM LMRMVLNESL KVLKMQGWEE RQIEELLELL MTEILHPSSQ APNGVKSHFI EIFLEELTKV GAEELTADQN LKFIDPFCRI AARTKDSLVL NNITRGIFET IVEQAPLAIE DLLNELDTQD EEVASDSDES SEGGERGDAL SQRSEKPPA GSICRAEPEA GEEQAGDDRD SGGPVLQFDY EAVANRLFEM ASRQSTPSQN RKRLYKVIRK LQDLGGIFP

image



Background

Alternative Names Ribosomal RNA processing protein 1 homolog A (Novel nuclear protein 1) (NNP-1) (Nucleolar protein Nop52) (RRP1-like protein)

Note

Product Name: Recombinant Human RRP1 Protein
Catalog #: PEH10306



For Research Use Only , Not for Diagnostic Use