


Product Name: Recombinant Human POU2F3 Protein
Catalog #: PEH10307



Summary

Name	PO2F3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UKI9
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YPSLSEVFGR KKKRSTIET NIRLTLEKRF QDNPKPSSEE ISMIAEQLSM EKEVVRVWFC NRRQKEKRIN CPVATPIKPP VYNSRLVSPS GSLGPLSVPP VHSTMPGTVT SSCSPGNNSR PSSPGSGLHA SSPTASQNNNS KAAVNSASSF

image



Background

Alternative Names	POU domain, class 2, transcription factor 3 (Octamer-binding protein 11) (Oct-11) (Octamer-binding transcription factor 11) (OTF-11) (Transcription factor PLA-1) (Transcription factor Skn-1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human POU2F3 Protein
Catalog #: PEH10307

