


Product Name: Recombinant Human ZNF521 Protein
Catalog #: PEH10309



Summary

| | |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name | ZN521_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q96K83 |
| Host | E.coli |
| Species | Human |
| Molecular Mass | 20 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | CKKRPEDGEE LEDEAVHSCD SCLQVFESLS DITEHKINQC QLTDGVDVED DPTCSWPASS PSSKDQTSPS HGEGCDFGEE EGGPGLPYPC QFCDKSFSRL SYLKHHEQSH SDKLPFKCTY CSRLFHKHRS RDRHIKLHTG DKKYHCSECD |

image



Background

| | |
|--------------------------|------------------------------------------------------------------------------------------------|
| Alternative Names | Zinc finger protein 521 (Early hematopoietic zinc finger protein) (LYST-interacting protein 3) |
|--------------------------|------------------------------------------------------------------------------------------------|

Note

For Research Use Only , Not for Diagnostic Use