


Product Name: Recombinant Human FOXN4 Protein
Catalog #: PEH10317



Summary

Name	FOXN4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96NZ1
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAIHRSMANP EELDKLISDR PESCRPGKP GEPEAPVLTH ATTVAVAHGC LAVSQLPPQP LMTLSLQSV LHHQVQPQAH LAPDSPAPAQ TPPLHALPDL SPSPLPHPAM GRAPVDFINI STDMNTEVDA LDPSIMDFAL QGNLWEEMKD EGFSLDTLGA FADSPLGCDL GASGLTPASG

image



Background

Alternative Names Forkhead box protein N4

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FOXN4 Protein
Catalog #: PEH10317

