



Product Name: Recombinant Human DBR1 Protein
Catalog #: PEH10322



Summary

Name	DBR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UK59
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ENNTLGSPAA SELLEHLKPT YWFS AHLHVK FAALMQHQAK DKGQTARATK FLALDKCLPH RDFLQILEIE HDPSAPDYLE YDIEWLTLR ATDDLINVTG RLWNMPENNG LHARWDYSAT EEGMKEVLEK LNHDLVKPCN FSVTAACYDP SKPQTQMQLI HRINPQTTEF CAQLGIIDIN VRLQKSKEEH HVCGEYEEQD DVESNDSGED QSEYNTDTSA LSSINPDEIM LDEEDEDSEI VSAHSGMNTP SVEPSDQASE FSASFSDVRI LPGSMIVSSD DTVDSTIDRE GKPGGTVESG NGEDLTKVPL KRLSDEHEPE QRKKIKRRNQ AIYAAVDDDD

image



Background

Alternative Names Lariat debranching enzyme (EC 3.1.-.-)

Product Name: Recombinant Human DBR1 Protein
Catalog #: PEH10322



Note

For Research Use Only , Not for Diagnostic Use