


Product Name: Recombinant Human SRSF11 Protein
Catalog #: PEH10324



Summary

Name	SRS11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q05519
Host	E.coli
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PRRHKKKKK DKDKERSRDE RERSTSKKKK SKDKEKDRER KSESDKDVKQ VTRDYDEEEQ GYDSEKEKKE EKKPIETGSP KTKECSVEKG TGDSLRESKV

image



Background

Alternative Names	Serine/arginine-rich splicing factor 11 (Arginine-rich 54 kDa nuclear protein) (p54) (Splicing factor, arginine/serine-rich 11)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use