


Product Name: Recombinant Human MACROD1 Protein
Catalog #: PEH10327



Summary

Name	MACD1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BQ69
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FLKGLSDKQR EEHYFCKDFV RLKKIPTWKE MAKGVAVKVE EPRYKKDKQL NEKISLLRSD ITKLEVDAIV NAANSSLLGG GGVDGCIHRA AGPLLTDECR TLQSCKTGKA KITGGYRLPA KYVIHTVGPI AYGEPSASQA AELRSCYLSS LDLLEHRLR SVAFPCISTG VFGYPCEAAA EIVLATLREW LEQHKDKVDR LIICVFLEKD EDIYRSRLPH

image



Background

Alternative Names	ADP-ribose glycohydrolase MACROD1 (MACRO domain-containing protein 1) (O-acetyl-ADP-ribose deacetylase MACROD1) (EC 3.1.1.106) (Protein LRP16) ([Protein ADP-ribosylaspartate] hydrolase MACROD1) (EC 3.2.2.-) ([Protein ADP-ribosylglutamate] hydrolase MACROD1) (EC 3.2.2.-)
--------------------------	--

Product Name: Recombinant Human MACROD1 Protein
Catalog #: PEH10327



Note

For Research Use Only , Not for Diagnostic Use