


Product Name: Recombinant Human SGO2 Protein
Catalog #: PEH10328



Summary

Name	SGO2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q562F6
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESNNLSADTP CATVLDKQHI SSPELNCNNE INGHTNETNT EMQRNKQDLP GLSSESAREP NAECMNQIED NDDFQLQKTV YDADMDLTAS EVSKIVTVST GIKKKSNKKT NEHGMKTFRK VKDSSSEKKR ERSKRQFKNS SDVDIGEKIE NRTERS DVLD GKRG AEDPGF IFNNEQLAQM NEQLAQVNEL KKMTLQTGFE QGDRENVLCN KKEKRITNEQ EETYSLSQSS GKFHQESKFD KGQNSLTCNK SKASRQTFVI HKLEKDNLLP NQKDKVTIYE NLDVTNEFHT ANLSTKDNGN LCDYGTHNIL DLKKYVTDIQ PSEQNESNIN KLRKKVNRKT EIISGMNHMY EDNDKDVVHG LKKGNFFFKT QEDKEPISEN IEVSKELQIP ALSTRDNENQ CDYRTQNVLG LQKQITNMYP VQQNESKVNK KLRQKVNRRKT EIISEVNHLD NDKSIEYTVK SHSLFLTQKD KEIIPGNLED PSEFETPALS

image



Background

Product Name: Recombinant Human SGO2 Protein
Catalog #: PEH10328



Alternative Names

Shugoshin 2 (Shugoshin-2) (Shugoshin-like 2) (Tripin)

Note

For Research Use Only , Not for Diagnostic Use